

# LADEBOY

Wheelchair loading systems



## On the way.

Without a vehicle the agility is set limits.

Most of all for wheelchair users. A reasonable integration to the daily routine is bound to the vehicle. And therefore to the necessity of loading and unloading of the wheelchair.

In a conventional way this is a certain challenge of physical strength and skills even for the accompanying person. Sometimes it is even an invincible barrier which only can be solved with the help of a third party. Therefore activities are limited needless or completely prevented.

The wheelchair loading system for the trunk. An idea which got something going which did not exist in this form up to now and therefore nothing comparable exists.



Find a way of yourself.



As a daily routine for wheelchair users who are able to walk some steps or passengers with handicap. Foldable wheelchairs with or without electric motor like e-fix and e-motion can be loaded easily.

The **LADEBOY** loads and unloads the wheelchair in the trunk of your vehicle without any effort. By pushing a button, without any strength, without bending down, without lifting. Just open the trunk, retract the **LADEBOY** with the remote control by pushing a button, move the wheelchair against it, lock the fastening handle and load automatically again.

The wheelchair is unloaded in reverse order. With the loading system the wheelchair is protected excellently from slipping, especially on rapid and curvy drives. A big relief.



No more lifting, let the **LADEBOY** do the loading.



Loading ...



The line of products of the LADEBOY convinces due to a modular construction kit. Many versions are possible. You have free choice of vehicles.

The different kind of physical handicaps request great flexibility to the loading system - individual solutions which are adapted almost perfectly to all needs. For the installation of this wheelchair loading system almost all types of vehicles are suitable. The number of seats in the vehicle can mostly be remained. An installation is possible on the right side as well as on the left side in the trunk.

**Useful.** ■ With quick fasteners the LADEBOY can be unloaded and loaded again at any time in a few minutes without any tools to use the trunk completely as needed.

The current consumption is made at the on-board socket or the cigarette lighter. A TÜV registration is not necessary. The LADEBOY is absolutely maintenance-free.



The current consumption is made at the cigarette lighter.



On shorter hatchback vehicles just a part of the seat bench is folded forward.

Which car can be used.



... and unloading again.

Intelligent customized installation kits take the fastening in the vehicle.

**Installation kit.** ■ The fastening of the LADEBOY is made with an installation kit, preferable at the available original screw attachment points of the vehicle. These guarantee a high security standard and a backfitting without any problems. The installation kit considers the particular type of vehicle; it is perfectly adjusted on it and contributes to the value preservation of the vehicle.



Installation on the left side possible.



**Version in VAN - optimum use of space.**

On vehicles with high trunk - VAN and Minivan - this advantage can be profited. The LADEBOY S convinces here with its compact design. The wheelchair can be loaded upright.



Weight classes are no barrier anymore.

**Strong performances.** ■ The LADEBOY is available as standard version with lifting power up to 20 kg and as maximum version with lifting power up to 50 kg. This means that wheelchairs with electric motor can also be loaded with the LADEBOY Maximum besides lightweight wheelchairs.

**Useful equipment. Split battery for e-fix.** ■ On wheelchairs with e-fix motor the battery can be removed for folding. It is easier in conjunction with the split battery. The wheelchair can be folded and loaded together with the battery. Details that inspire. Further equipment is the additional control unit. It makes sure that you have even more ease of use while loading heavy wheelchairs.

**Quality.** ■ A fine aid is recognized by the function, the design and the high-quality components. But only with a continuous further development a completely perfected product is developed. In the last 10 years the LADEBOY has been refined. Snow and minus temperatures in Norway, sand and heat in Israel. The LADEBOY works reliable. Especially for wheelchair users safety, maintenance-freedom and stable value are of particular importance.



You have the choice.



## Expeditions.

Go shopping. Visit friends. Create freedom. Go out in all weathers.

A vehicle with a loading system is a worthwhile investment for a self-determined life.

In doing so enjoy the ease of loading, be happy about clean clothes and spare your body. Make a healthy decision.





#### Quality from responsibility.

High-quality materials for reduced friction.

Very little current consumption: 0,04 Ah.

Power cable pluggable at the on-board socket or cigarette lighter.

Construction kit enables solutions for almost every trunk design, as well in a cabriolet.

Compact construction: width 30cm, height 11cm, length from 82cm up, depending on the vehicle.

Usually reusable when changing vehicle or backfitting.

Reliable technique; a thousand times proven linear technology.

New from model 2006 up:

- "Pro Secure Technologie" for more safety for the passengers.
- Emergency release.
- Fast loading speed.
- Maintenance-free.



#### ■ LADEBOY Standard

Suitable for folded wheelchairs up to 20 kg weight of wheelchair.

#### ■ LADEBOY Maximum

Suitable for folded wheelchairs and wheelchairs with electric motor (e-fix, e-motion) up to max. 50 kg weight of wheelchair.

To the specifications in this catalogue. After the deadline of this brochure, 01.03.2006 modifications can arise at the product. Design or dimensional modifications, colour deviations on the part of the manufacturer remain reserved during the delivery time in case the modifications or deviations are reasonable for the customer in consideration of the vendor's interests. Depending on the type of the vehicle certain measurements of the wheelchair like length, width and height are necessary due to the vehicle's dimensions, seats can be omitted. On request we send you the particular measurements under specifications of the type of your vehicle. LADEBOY is a registered trademark. Usage of the trade name or extracts out of this brochure only with written approval. Registered at the German Patent Office, Rausch Technik GmbH

Wheelchair loading systems

RAUSCH TECHNIK GMBH

Ind.Geb. Lehenmorgen

Hölzlestraße 27

D-72336 Balingen

Phone +49(0)7433-808-1

Fax +49(0)7433-808-2

info@rausch-technik.de

www.LADEBOY.de

TECHNIK  
**RAUSCH**